


HEINE QUALITY
MADE IN GERMANY



∴ HEINE ML 4 LED HeadLight with DV1

Digital Video Camera
for high-resolution images.

HD

HIGH-RESOLUTION
DIGITAL IMAGES

**LIGHT
WEIGHT**

LIGHTWEIGHT



OPTIMAL
INTEGRATION



HEINE® ML 4 LED HEADLIGHT WITH DV 1. DIGITAL VIDEO CAMERA FOR HIGH-RESOLUTION IMAGES.

High-resolution digital images of every examination. We make it happen with the new HEINE DV 1 Digital Video Camera and its smooth integration in the coaxial illumination of the ML 4 LED Headlight. Precision HEINE Optics guarantees that even the smallest details can be seen on the camera image. The compact camera is adjusted to the field of illumination in just a few seconds. Due to the light weight of the camera, the ML 4 LED HeadLight provides the usual comfortable and secure fit even in combination with the DV 1. A standard USB 2.0 interface ensures reliable and uncomplicated data transfer.

HD

HIGH-RESOLUTION DIGITAL IMAGES

High camera resolution of 1,280 x 960 Pixel is perfect for sharing examination images or examination techniques with a larger audience.

LIGHT WEIGHT

LIGHTWEIGHT

The compact camera is so well integrated that its light weight is barely noticeable in combination with the ML 4 LED HeadLight.



OPTIMAL INTEGRATION

Ideal concept for a perfectly matched system: ML 4 LED HeadLight, DV 1 Digital Video Camera and optional HEINE Binocular Loupe.



ML 4 LED HeadLight with DV 1

ML 4 LED HeadLight with DV 1 Digital Video Camera and mPack UNPLUGGED with plug-in transformer [01]

J-008.31.470

ML 4 LED HeadLight with DV 1 Digital Video Camera and mPack with plug-in transformer [02]

J-008.31.471

